

RNF148-ring finger protein 148 Gene

Reference: PTXBC029264

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RING finger protein 148		
For research use only.			

GENE DESCRIPTION	ring finger protein 148
GENE ID	378925
ORIGIN SPECIES	Human
ACCESSIONS	BC029264
VECTOR	pDONR223
NCBI SYMBOL	RNF148
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcttcttagaattacccttcgacgcatgttctgttcatctggacttttgaggcttagtatcttctactacttgccttgcactaaacggaagccatttggacagctcacctgaataacatttcaggttggaaatgagatcacatcggaattaggagagtgagtggttcggga
atcattctctctggaaagggtgtctgtgtgtgtgacacttctgaaggatggaatcagaatgcctgtcatccttgcaccaatttcagcaggcccaaacagcgagactcatggctggccctcatcgaacgtggaggctgtacttttacataaaatcaacgtggcagcagagaaggagcaa
atgggtgatcatctacaactatcaagggtacggcgagtaaaagtatttccatgtctcaccaggggacgggaaaatagtcgcggtgatgataagcaacctgaaaggcatggaaatttgcactcgattcagaaggagtgatgtgacagtcattgaagtggggagaatgcacatgcagt
gggtgagccattacatcatgtatctatttaccttctgctgctccacaattgcctacttttagattgcgtctggagacttacacctagagtgcccaattcttccaccaggaggcgaagtcaataaagacagatgtgaagaaagctattgaccagcttcaactcgcagttctcaagaaggga
tgagggaattagacctaaatgaagacaactgtgtgttgccttgacacataaaccccaagatgtagtagcattttaacttgcaaacatttttccataaggcatgcatgacccctggcttttagcccataggacatgtccatgtgcaagtgtgacatcctgaaacttaa