

ATG10-ATG10 autophagy related 10 homolog (S. cerevisiae) Gene

Reference: PTXBC029268

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATG10 autophagy related 10 homolog; ubiquitin-like-conjugating enzyme ATG10; APG10; APG10L; pp12616; APG10 autophagy 10-like; autophagy-related protein 10; autophagy related 10		

For research use only.

GENE DESCRIPTION	ATG10 autophagy related 10 homolog (S. cerevisiae)
GENE ID	83734
ORIGIN SPECIES	Human
ACCESSIONS	BC029268
VECTOR	pDONR223
NCBI SYMBOL	ATG10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagaagatgagttcattggagaaaaacattccaacgttattgtgcagaattcattaacattcacacagataggtgatagttgggaatggagaccatcaaggactgttctgatggctacatgtgcaaaatacattcaaatgaagaatgggtctgtgatgtcacatctaggagcatcta
cccatggacagacatgtctccatggaggaggtttcgagctaccccttggatgattgtgaagtgtgaaactgcagcagcgtccgaagtgtatgaatgatgatcatgtcttattctgtagctaccaagtgctgtacttactttagggcaagcttttagatgggagaccttaactctga
aggacatatgggaaggagttcatgagtgctataagatgcgactgtacagggaccatgggacactattacgcaacagggaacatccaatactgggcaaccctttttgtacttcatcctgcaagacgaatgaattcatgactcgttattaagaatttcagaaaaataagaatgtcaact
atatcacatcatggctgagcattgtagggccagttgtgggtgaatctacctctgagttatgccaaagcaacgtctcaggatgaacgaaatgtcccttaa