

RBBP6-retinoblastoma binding protein 6 Gene

Reference: PTXBC029352

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RBBP6; MY038; P2P-R; PACT; RBQ-1; SNAMA; RB-binding Q-protein 1; p53-associated cellular protein of testis; proliferation potential-related protein; protein P2P-R; retinoblastoma binding protein 6; retinoblastoma-binding Q prot		

For research use only.

GENE DESCRIPTION	retinoblastoma binding protein 6
GENE ID	5930
ORIGIN SPECIES	Human
ACCESSIONS	BC029352
VECTOR	pDONR223
NCBI SYMBOL	RBBP6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctgtgtgcattataaatttctctaaactcaactatgataccgtcaccttggggctccacatccctctgcgactaaagaagcagattatggggagagagaagctgaaagctgccgactgcgacctgcagatcaccaatgcgcagacgaaagaagaatactgatgataatgc
tctgattcctaagaattcttctgttaattgttagaagaattcctattggaggtgttaaatctacaagaagacatagtttataatgcgaactgaaccagcagatggcaactacaaaagcagtatgtaaaacacaatctcacactttttacacattgcttttacccttataa