

## RBP5-retinol binding protein 5, cellular Gene

Reference: PTXBC029355

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRBP-III; CRBP3; CRBP3III; HRBPiso; retinol-binding protein 5; cellular retinol-binding protein III; retinol binding protein 5, cellular; retinol binding protein 5		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	retinol binding protein 5, cellular
<b>GENE ID</b>	83758
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC029355
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	RBP5
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atgctcccaacctcactggctactaccgcttctcgcagaagaacatggaggactacctgcaagccctaaacatcagcttggctgtcgcgaagatcgcgtgctgctgaagccggacaaggagatgaacaccaggccaaccatgacggtgaggacgctcagcacctccgaaa  
ctacactgtgcagttgatgtggagtgagttgaggaggacctcaggagcctggacggcggacgaaaatgccagaccatagtaacctgggaggaggagcacctggtgtgtgtgcagaaagggaggtcccaaccgggctggagacactggctggaggagagatgctgtatctg  
aactgactgcaagggatgagctgtgagcaggtctcaggaaggtcagatag