

TPSB2-tryptase beta 2 Gene

Reference: PTXBC029356

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TPS2; tryptaseB; tryptaseC; tryptase beta-2; mast cell tryptase beta II; mast cell tryptase beta III; testicular tissue protein Li 163; tryptase beta 2 (gene/pseudogene)		

For research use only.

GENE DESCRIPTION	tryptase beta 2
GENE ID	64499
ORIGIN SPECIES	Human
ACCESSIONS	BC029356
VECTOR	pDONR223
NCBI SYMBOL	TPSB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgaatctgtctgtggcgtgcccgtctgtgagcggcgccctacggcgccctgccccagggccctgcagcagtggtggcctgttgggggtcaggaggccccaggagcaagtggccctggcaggtgagcctgagagtcacggccatactggatgcacttctg
cgggggtcctcatccacccagtggtgtgacccgacgacgtggtgacggacgtcaaggtatctggccgctcagggtgcaactgcgggagcagcactctactaccaggaccagctgtgctgggtcagcaggtatcgtgcacccacagttctacacggccagat
cggagcggacatcgccctgtgtgagctggaggagccgtggaaggtctccagccacgtccacacgggtcacctgccccctgctcagagacctccccccggggatgccgtgtgtgtgctggtgctggggcgtggtgacatgatgagcgcctcccaccgcatctctgaagca
ggtgaaggtccccataatgaaaaccacatttgtgacgaaataccacgttggcgctacacgggagacgacgtccgcatctgctgacgacatgtgtgtgcccgggaacccggaggaggtactatgccaggggcactccggaggccctgtgtgcaaggtgaatggcacctg
gctgcaggcggcgtgtgacgtgtgggcgagggctgtgcccagccaaccggcctgacatctacacccgtgtacactactactgtgactggatccaccactatgtcccaaaagccgtga