



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SUV39H2-suppressor of variegation 3-9 homolog 2 (Drosophila) Gene

Reference: PTXBC029360

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	histone-lysine N-methyltransferase SUV39H2; KMT1B; H3-K9-HMTase 2; histone H3-K9 methyltransferase 2; lysine N-methyltransferase 1B; su(var)3-9 homolog 2; suppressor of variegation 3-9 homolog 2		

For research use only.

GENE DESCRIPTION	suppressor of variegation 3-9 homolog 2 (Drosophila)
GENE ID	79723
ORIGIN SPECIES	Human
ACCESSIONS	BC029360
VECTOR	pDONR223
NCBI SYMBOL	SUV39H2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]