

OCEL1-occludin/ELL domain containing 1 Gene

Reference: PTXBC029361

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FWP009; S863-9; occludin/ELL domain-containing protein 1; occludin/ELL domain containing 1		

For research use only.

GENE DESCRIPTION	occludin/ELL domain containing 1
GENE ID	79629
ORIGIN SPECIES	Human
ACCESSIONS	BC029361
VECTOR	pDONR223
NCBI SYMBOL	OCEL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacaacccggacggaagtgcctctccgacagcagatccaggctcggagctccagacgtgggacagcgcccgagaccaccccgccgcgccgggacacgacgcccccgaggaccccgccatcagccggaaaccctgagctgcttctccggaggccgatgcc
cacccgggagcccccaagactcgcggctcccgggggcactgcatactcaccgctgggctgggccccgctgcagggactggcgccccgaggcctcaaaaccagcgcccccgccctccgtgccagccccagccgggacccccacaaggcaagaccaagaagattgtgt
ttgaggatgagttgctctccagccctctggcgccaagaagcctattggagccatccctaaggggcataagcctaggccccaccagtgcccgactatgagcttaagtaccgcccagtgagcagtgagagggaacggagccgctatgtcgcagtggtccaggaccagtaaggaga
gttcttgagctccagcacgaggtgggtgtgacagggcaaaagtcaggcagctggaggccctgctgagctcctgccccaccacaaagcagaaggaggcccaagttgcagccgggtttggaggaggtttgagatgaagcgaatggatcctggtctcctggacaagcaggctcg
ctgccactactgaagggtaaactgaggcatctcaagactcagatccagaaattcagtgaccaaggagacagcagggctcgtgtacttctaa