

GBA3-glucosidase, beta, acid 3 (cytosolic) Gene

Reference: PTXBC029362

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CBG; CBGL1; GLUC; KLRP; cytosolic beta-glucosidase; cytosolic beta-glucosidase-like protein 1; glucosidase, beta, acid 3 (cytosolic); klotho-related protein; glucosylceramidase beta 3 (gene/pseudogene)		

For research use only.

GENE DESCRIPTION	glucosidase, beta, acid 3 (cytosolic)
GENE ID	57733
ORIGIN SPECIES	Human
ACCESSIONS	BC029362
VECTOR	pDONR223
NCBI SYMBOL	GBA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttccctgcaggatttgatggcgccagccactgcagcttatcaagtagaaggagctgggatgcagatggaaaaggccctgtgtctgggacacatttactcatcagggaggagagagagtttcaagaaccagactggcgatgtagctgtggcagctacactctgtgggagg
aagattgaaatgatcaaacagcttgatgactcattaccgtctctcttctcgtgacgtctgttacctgatgggacgacaggtttcatcaaccagaaagctatccaacttgataaaagtcaatcttcaagtataattgtgcatggtctctcttgataactttgagtggaaccaggatcacagcag
ccggttgggtctctccacgttgatttgaagaccagctagaccccgagtccttacatcggccaaggaatatgccaagatcatccgaacaatggccttgaagcacatctgtag