

GNGT1-guanine nucleotide binding protein (G protein), gamma transducing activity polypeptide 1 Gene

Reference: PTXBC029367

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GNG1; guanine nucleotide-binding protein G(T) subunit gamma-T1; guanine nucleotide binding protein (G protein), gamma transducing activity polypeptide 1; transducin gamma chain; G protein subunit gamma transducin 1		

For research use only.

GENE DESCRIPTION	guanine nucleotide binding protein (G protein), gamma transducing activity polypeptide 1
GENE ID	2792
ORIGIN SPECIES	Human
ACCESSIONS	BC029367
VECTOR	pDONR223
NCBI SYMBOL	GNGT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagtaatcaatattgaggacgtgacagaaaaggacaaattgaagatggaagttgaccagctcaagaaagaagtgacactggaaagaatgctagttccaaatgttgaagaagtaagagattacgttgaagaacgatctggcaggatccactggtaaagggcatccagaggac
aaaaatccctcaaggagctcaagaggagctgtgtgattcataa