

GGA1-golgi associated, gamma adaptin ear containing, ARF binding protein 1 Gene

Reference: PTXBC029388

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ADP-ribosylation factor-binding protein GGA1; ADP-ribosylation factor binding protein 1; gamma-adaptin related protein 1; golgi-localized, gamma ear-containing, ARF-binding protein 1; golgi associated, gamma adaptin ear containing, ARF binding protein 1		

For research use only.

GENE DESCRIPTION	golgi associated, gamma adaptin ear containing, ARF binding protein 1
GENE ID	26088
ORIGIN SPECIES	Human
ACCESSIONS	BC029388
VECTOR	pDONR223
NCBI SYMBOL	GGA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcccgcatggagccggagactctggagcgcgaaatcaatagagccacgaacccctgaacaaggagctcgactgggccagatcaacggctctgcgagcagctcaacgaggactttgaggggctcactcgcacccggctgtggtgccacaagatccagtccccc
aggagtgggagggcatccaggcccttgacggtgagaaggggagagggccaccatccgtcccccgcctgtgacgacaccaaaggagggccaagactga