

PSMA3-proteasome (prosome, macropain) subunit, alpha type, 3 Gene

Reference: PTXBC029402

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HC8; PSC3; proteasome subunit alpha type-3; macropain subunit C8; multicatalytic endopeptidase complex subunit C8; proteasome (prosome, macropain) subunit, alpha type, 3; proteasome component C8; proteasome subunit C8; testicular secretory protein Li 43;		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) subunit, alpha type, 3
GENE ID	5684
ORIGIN SPECIES	Human
ACCESSIONS	BC029402
VECTOR	pDONR223
NCBI SYMBOL	PSMA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctcaatcgccactgggtatgacctgtcagcctctacattctctctgacggaagagttttcaagttgaatatgctatgaaggctgtgaaaaatagtagtacagctattggaatcagatgcaaatggtgtgttctttgggttagaaaaatagtcctttctaaactttatgaagaaggttcc
aacaaaagactttttaaattgtgatcgccatgttggaatggcagtagcaggtttgttgccagatgctcttcttagcggacatggcaagagaagaagcttcaacttcagatcctaaactttggctacaacattccactaaaacatcttgcagacagatggccatgtatgtcatgcataatacactta
cagtgtctgttagacttttggctgcagttcatgttaggtcttactagtgtaagacggtgacgaactctacatgattgacccatcaggtgtttcctacaggttattggggctgtgcccacgcaaaagccagacaagctgcaaaagacggaaatagagaagcttcagatgaagaagaatgacctgcc
gtgatatcgttaaagaagttgcaaaaataattacatagtagacgaagttaaggataaagctttgaactagaactcagctgggtgtgtgaattactaatgaagacatgaaattgttccaaaagatataagagaagaagcagagaatatgctaaggaatctctgaaggagaagatgaa
tcagatgatataatgtaa