

H3F3A-H3 histone, family 3A Gene

Reference: PTXBC029405

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H3.3A; H3F3; histone H3.3; H3 histone, family 3A; H3 histone family member 3A		
For research use only.			

GENE DESCRIPTION	H3 histone, family 3A
GENE ID	3020
ORIGIN SPECIES	Human
ACCESSIONS	BC029405
VECTOR	pDONR223
NCBI SYMBOL	H3F3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcgtacaaagcagactgccgcaaatcgaccgtggtgtaagcaccagggaagcaactggctacaaaagccgctcgcaagagtgccctctactggaggggtgaagaaacctcatcgttacaggcctggtactgtggcgctccgtgaaattagacgttatcagaagtcactga
acttctgattcgcaaaacttccctccagcgtctgtgctgagaaattgctcaggactttaaaccagatctgcgttcagagcgcagctatcgggtcttgcaggaggcaagtgaggcctatctggttggcctttttgaagacaccaacctgtgtgctatccatgccaaacgtgaacaattatgcc
aaaagacatccagctagcagccgcatactggagaacgtgcttaa