

GPX8-glutathione peroxidase 8 (putative) Gene

Reference: PTXBC029424

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EPLA847; GPx-8; GSHPx-8; UNQ847; glutathione peroxidase 8 (putative)		
For research use only.			

GENE DESCRIPTION	glutathione peroxidase 8 (putative)
GENE ID	493869
ORIGIN SPECIES	Human
ACCESSIONS	BC029424
VECTOR	pDONR223
NCBI SYMBOL	GPX8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcctcttgcagcttaccgctaaaatgtccgggccagagcaaaaggtatttgcagtttctgtctatagttctatgcacagtaacgctatttcttacaactaaaatcctcaaacctaaaatcaacagcttttatgcctttgaagtgaaggatgcaaaaggagaactgtttctctggaaa
agtataaaggcaaaagtttactagttgtaaacgtggccagtgactgccaactcagacagagaaattacttagggctgaagggaactgcacaaagagtttgaccatccacttcagcgtgttggtttccctgcaatcagtttgagaaatcggagcccccccaagcaaggaagtagaatcttt
tgcaagaaaaactacggagtaactttcccatctccacaagattagattctgaaggagaacctgcatfttagatttctgttgattctcaagaaggaaaccaagggtggaattttggaagtatctgtcaacctgaggggtcaagttgtgaagttctggaaggccagaggagccatt
gaagtcacaggccctgacatagcagctctggttagacaagtgatcataaaaaagaagaggatctatga