

CCDC23-coiled-coil domain containing 23 Gene

Reference: PTXBC029427

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCDC23; small vasohibin-binding protein; coiled-coil domain containing 23; coiled-coil domain-containing protein 23; small vasohibin binding protein		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 23
GENE ID	374969
ORIGIN SPECIES	Human
ACCESSIONS	BC029427
VECTOR	pDONR223
NCBI SYMBOL	CCDC23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatccacgtgcacgtaaagaaaaaaccaagtttaagaatctgtcagcagagttgagaaggccaacagaaatcagcccagcaggagctgaagcagagacaaagagcagagatctatgccctcaacagagtcacagaaactggagcagcagcagtttgatgagttctgtaaac
agatgcagcctctggagaatga