

APITD1-apoptosis-inducing, TAF9-like domain 1 Gene

Reference: PTXBC029430

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APITD1; CENP-S; FAAP16; MHF1; centromere protein S; FANCM associated histone fold protein 1; FANCM-interacting histone fold protein 1; Fanconi anemia-associated polypeptide of 16 kDa; apoptosis-inducing TAF9-like domain-containing protein 1; apoptosis-ind		

For research use only.

GENE DESCRIPTION	apoptosis-inducing, TAF9-like domain 1
GENE ID	378708
ORIGIN SPECIES	Human
ACCESSIONS	BC029430
VECTOR	pDONR223
NCBI SYMBOL	APITD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggaggaggcggagaccgaggagcagcagcattcttaccacagaggctaaggcagcaggttactatctgtgggtgtcttgcgaggaagtgcatggacaagagatgcagttcagcaaacagaccattgcggccatttcggagctgactttccagcagtgtaaaatttgcacaaagaccttgaaatgtttgcaagacatgcgaaaagaaccacaaataaacctgaagatgtgaagctcttagccaggaggagtaattcactgctaaaatacatcacagacaaaagtgaagagattgctcagattaacctagaacgaaagcacagaaaaagaagtcagaggatgaagcaaaaattcaaggcagccagcagagcgtggagtgtggaaagtgaataaa