

## ZMYM6-zinc finger, MYM-type 6 Gene

Reference: PTXBC029439

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Buster2; MYM; ZBED7; ZNF198L4; ZNF258; zinc finger MYM-type protein 6; transposon-derived Buster2 transposase-like protein; zinc finger protein 258; zinc finger, BED-type containing 7; zinc finger, MYM-type 6; zinc finger MYM-type containing 6		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	zinc finger, MYM-type 6
<b>GENE ID</b>	9204
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC029439
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	ZMYM6
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atggcgtctctgtcttgcgtgctgctgctgctggcgcctctctcgcgctcgtggggtggccaagctcggaggagatcggctccagttcggagcggatgaatgccctgttcgtgcagtttgctgaggtgtcccgtgaaggtattggctaccagccagatcccctgaactacca  
aatagctgtggccttctggaactgctgctgggtgctgctgctatggccaccgatcgaagagatcagtaactgttcttgattctcctatgatggggctatcttcacctggcagctcgaagagatcactaagcactgtatcccagccattgtcgtcgggttctgctgctg  
ctgaatgtcggccagctcttagcccagactaagaaggtgctcagaccactaggaagaagactctaagtacattcaaggaatcctggaagtag