



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

EHP1-EH domain binding protein 1 Gene

Reference: PTXBC029477

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HPC12; NACSIN; EH domain-binding protein 1; NPF calponin-like protein; testis tissue sperm-binding protein Li 50e; EH domain binding protein 1		

For research use only.

GENE DESCRIPTION	EH domain binding protein 1
GENE ID	23301
ORIGIN SPECIES	Human
ACCESSIONS	BC029477
VECTOR	pDONR223
NCBI SYMBOL	EHBP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttcagtttgggaagagactgcagcgtgtgggaaacatgcatccaagtgtccagtttggccctcaccaggagctcatggttgagtgtacgaagaataatgcaaccagataaactgtgggtagtttggaccagaagaagccgaaggaaagtctctaaggcacatagctggcaacctg
gaaataaaaaatccctatcgtgctgtttgtgtgctgctgttctcgaacaattgaaatcactgttaacacttttaagatccctatgcggaagaattgaagacaaagatggacatttgcatagaatagttaagctaa