

CLRN3-clarin 3 Gene

Reference: PTXBC029478

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM12; USH3AL1; clarin-3; transmembrane protein 12; usher syndrome type-3A-like protein 1; clarin 3		

For research use only.

GENE DESCRIPTION	clarin 3
GENE ID	119467
ORIGIN SPECIES	Human
ACCESSIONS	BC029478
VECTOR	pDONR223
NCBI SYMBOL	CLRN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctaccacaaagaagacattgatgttctatcaagctttaccagccttgggtccttcattgtaattgtctatcttgggacacaagcatggatcaccagtagtgctgttagagactctgttcaaatgggagcatttcatcacttacggacttttctgtgggagagtagtgaagaatt
gagtcacggacttgcagaaccaagaaaaagttgacgttttagagatactgaataattctcccaaaaactctgcatcgggtgactatctgttctgtgctgagtttgatcacgtcgtgctgagctctgggtttacctctacaacagcatcagcaaccttaccagacattcctgggccc
acgggggtgtacacctggaacgggctcgtgcatcctctgttttggacatgatactgttggcgaacacgcagtcaccaactctcgaagagttgttccaaatgctttaccgggcaaccaccagtaaggaacgacccacagttacggatactgttctggctcactgtctgcat
tcttctaataatagtcactgttaacatcatcattttctaccagaaggccagataccagcgggaagcaggagcagagaaagccaatggaatatgctccaaggagacggaattttattctga