

TPRKB-TP53RK binding protein Gene

Reference: PTXBC029492

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EKC/KEOPS complex subunit TPRKB; CGI-121; PRPK (p53-related protein kinase)-binding protein; TP53RK binding protein		

For research use only.

GENE DESCRIPTION	TP53RK binding protein
GENE ID	51002
ORIGIN SPECIES	Human
ACCESSIONS	BC029492
VECTOR	pDONR223
NCBI SYMBOL	TPRKB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagttacacatcagctggacatttccgaatgcagggtaccctctgtatttaagatgtaaaaaatcgaggagacttgagaagaaggccatggaaggcaccatcgatggatcactgataaatcctacagtattgttgatcatttcagatactgtggcagcaacaaagcagtt
cacctctacaaactgggaaaatgaagacaagaactctatctactgaaattatttcaaccttccccaataacaatttcagaggcttgaaaaaattggtatctcagcaaatgacacttcaattcttaattgtttacattgaaggaggagaaaaaataatacaagaatacctaataatctcaag
tagaaggatcatcaggtttctctgaaaaatcttctgaaataatgaattatcagaagtcacaaagatatataaactctcttcacaagaagaagattgggacattattggatgctatcattgtagaatgtcaacaaagatgtttatga