

PTPRS-protein tyrosine phosphatase, receptor type, S Gene

Reference: PTXBC029496

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PTPSIGMA; R-PTP-S; R-PTP-sigma; receptor-type tyrosine-protein phosphatase S; protein tyrosine phosphatase PTPsigma; receptor-type tyrosine-protein phosphatase sigma; protein tyrosine phosphatase, receptor type S		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase, receptor type, S
GENE ID	5802
ORIGIN SPECIES	Human
ACCESSIONS	BC029496
VECTOR	pDONR223
NCBI SYMBOL	PTPRS
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggccacgtggggccctggcatggtgtctgtgtgtgtccatgggcccctctgtgtctctgtctgtgtggagctgtgcagcagaagagccccaggtttatcaagaaccaaggaccagatcggtgtgtcgggggtgtggcctcttctgtgtcaggccacgggtgacccc
aagccacgagtgacctggacaagaagggaaggaaggtcaactctcagcgtttgagacgattgagtttgatgagagtgccgggagtgctgaggatccagccgtgaggacaccggtgatgaaacgtgtacgagtggtgtggccagaactcggttgggagatcacagtcct
gccaagcttactgtctccgaggacctga