

NPY-neuropeptide Y Gene

Reference: PTXBC029497

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PYY4; pro-neuropeptide Y; prepro-neuropeptide Y		

For research use only.

GENE DESCRIPTION	neuropeptide Y
GENE ID	4852
ORIGIN SPECIES	Human
ACCESSIONS	BC029497
VECTOR	pDONR223
NCBI SYMBOL	NPY
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctaggtaacaagcgactgggctgtccggactgacctcgcctgtccctgtcgtgtgcctgggtgcgtggcgaggcgtaccctccaagccggacaacccggcgaggacgcaccagcggaggacatggccagatactactcagcgctgcgacactacatcaacctcatc
accaggcagagatatggaaaacgatccagccagagacactgatttcagacctcttgatgagagaaagcacagaaaatgttccagaactcggcttgaagaccctgcaatgtgtgtga