

TEX12-testis expressed 12 Gene

Reference: PTXBC029506

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	testis-expressed sequence 12 protein; testis expressed sequence 12; testis-expressed sequence 12 protein variant 1; testis-expressed sequence 12 protein variant 2; testis-expressed sequence 12 protein variant 3; testis expressed 12		

For research use only.

GENE DESCRIPTION	testis expressed 12
GENE ID	56158
ORIGIN SPECIES	Human
ACCESSIONS	BC029506
VECTOR	pDONR223
NCBI SYMBOL	TEX12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggcaaatcaccctgtaaagcctgataataggaattgcaagaggccaagagaattggagtcctccagtgccagatagtcacagctgtcctctctggaaaacagattcatcttctctgaaattccggactattttataagatgaagccttgagaaaagattaaatgatgtgagcaagg
aaattaatctaattgtgtctacatgcaaaagcttttaagtgagagagcagcagtagatgcattcttaccattgatgatagatgaactcttcaagaagccaatgctattgaaaactttctaatacaaaaaagagagttctcgacagaggtttacagtgattgcaaacacattacacagataa