

DNAJB8-DnaJ (Hsp40) homolog, subfamily B, member 8 Gene

Reference: PTXBC029521

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT156; DJ6; dnaJ homolog subfamily B member 8; DnaJ (Hsp40) homolog, subfamily B, member 8; testicular tissue protein Li 56; DnaJ heat shock protein family (Hsp40) member B8		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily B, member 8
GENE ID	165721
ORIGIN SPECIES	Human
ACCESSIONS	BC029521
VECTOR	pDONR223
NCBI SYMBOL	DNAJB8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctaactactacgaagtgcctggcgctcaggccagccttccccggaggacatcaagaagcctaccgcaagctggccctcgttggcaccccgaagaacctgacaataaggaggaggcggagaagaagttcaagctggtctgaggccctatgaggtctgtgactcaa
gaaacgtccctgtatgacctgctggcctgtgacagctggcggtgctggcggtggccagcagccctaccacagcccttcgacacggcctacaccttcgtaacctgaggacatcttccggagtttttgggtggcctggacctttctctttagttctgggacagccattcaatag
tgacctggtggcggtggccatggcctgagggggctctcggcaggctttggagaattccggccttcattgaggccttctcatcttcaacatgctggcctgcagcggggcgagccacacccttctcatcctcctctggggctcagttctggcagctcgggttcaagtcg
gtgatgtctccacgagatgatcaatggccacaagtcaccaccaagcgcctcgtggagaacggcgaggcgcgtggaggtggaggagacggcgagctcaagtcggtgactgtgaacggcaaggagcagctcaatggatggacagcaagtag