

PFN4-profilin family, member 4 Gene

Reference: PTXBC029523

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	profilin-4; profilin IV; profilin family member 4		

For research use only.

GENE DESCRIPTION	profilin family, member 4
GENE ID	375189
ORIGIN SPECIES	Human
ACCESSIONS	BC029523
VECTOR	pDONR223
NCBI SYMBOL	PFN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccattgcagacttattgtagacacctctgggaaccaagcatgtggacagtcagccctcatcaaaatccaggagcggagcctgtgtgtagcatcaccagtttcaatgtaacgccagtgatgccgaacctgtgaatggattccaagaacctttgcaagcccgaaga
gaaggactgtattcaagggaagattacagatgttcggcagatgaattctctttagcctcaaaatgagaacctgtgtgtgtcgtgaagaccatctgtatctctgtagcaacttacctgaggcatgtatcctagcatctgttgaagccacagagagcctgggagacta
cctaagaaaaaaggaaagttaa