

## MAGEB18-melanoma antigen family B, 18 Gene

Reference: PTXBC029525

|                |                                                                                                                                       |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                   | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                  | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | melanoma-associated antigen B18; MAGE-B18 antigen; melanoma antigen family B, 18; melanoma antigen family B18; MAGE family member B18 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | melanoma antigen family B, 18                                                                     |
| <b>GENE ID</b>          | 286514                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC029525                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | MAGEB18                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgcctcaggtcagaagagtaagctccgtgcccgtgagaaacgccaccaggctcgttgagaaatcaggtatctggagctacgcaggccactgtggcagaaggagagtcaccctcctctgctatcttctttgtgacagacccagaattgcctgctgctgagacacatagcatcc  
ctgaagcgttcaggaggcccatccaccaccaatgctattgcacctgtttcatgcagttcaaatgaagggtgccagcagccaagatgagaaaagtctaggttcctcaagggaagctgagggctggaagaagatcctttaacaagaagtagtgcgctggtgctatttctgctcagaagt  
atgaaacgaagagccaattacaaggaggagataataaagtgtttatcagggaaggataaggtgctacttcaatgagatcctcaagagagcctctgagcacatggagctggcacttgggtgttgaattgaaggaaagtgatccatcaggcactactatgctttttcagcaaatagacctcacc  
tatgatgaaacaaccagtgatgaagaaaaattcccaagactggcctcctgatgattgcactgggtgtgatctttctgaatggcaaccgtgcccagaagagcagctctgggaattatgaatatgaggggtgtatgctgatagggaagcacttctctatggggatcccaggaaaggtcatga  
ccaaagatttggcagctaaagtacctggagtagcagcaagtgcacaacagtgatcctccacgtatgaattctgtgggtccaaagagctcacgctgaactagcaagatgaaagtcctggagttgtagccaagatacatgataccgtccctagtccttcccatcctgctatgaagagg  
ctttgaggatgagggaacagagaaccaagccagagctgcagccagggtctatctgctccatggcaaatgcacgttccagaaccacgtctagcagcttctcccatgctaagtga