

NUDT9P1-nudix (nucleoside diphosphate linked moiety X)-type motif 9 pseudogene 1 Gene

Reference: PTXBC029544

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	nudix (nucleoside diphosphate linked moiety X)-type motif 9 pseudogene 1
GENE ID	119369
ORIGIN SPECIES	Human
ACCESSIONS	BC029544
VECTOR	pDONR223
NCBI SYMBOL	NUDT9P1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggatccaggagagaaaattagtctacattgaagagagaatttggcagggaagccatgaactctttacagaaatctaggaaagaatgcaggagttggagagacaattgcacaaactcctcagtcaggaaatttgggtctataaaggctacgtggatgattctagaacacagat
aactgctggatagagacagaggctgtgaactaccgtgatgaaataatggaccatcttctctggaagctggcggatgccaagaaatgaaatgggtggacattaatgataaacttgagctttatgccagcctctcccaatttattcaacttggctgagaaacgaggtgccactga