

CLEC1B-C-type lectin domain family 1, member B Gene

Reference: PTXBC029554

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1810061I13Rik; CLEC2; CLEC2B; PRO1384; QDED721; C-type lectin domain family 1 member B; C-type lectin-like receptor 2; CLEC-2		

For research use only.

GENE DESCRIPTION	C-type lectin domain family 1, member B
GENE ID	51266
ORIGIN SPECIES	Human
ACCESSIONS	BC029554
VECTOR	pDONR223
NCBI SYMBOL	CLEC1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggatgaagatgatacatcaccttaafatataaaactcggaaaccagctctcgtctccgtctcatgcagcgcattacctaagaatgagaatgaaaatgcacaggaactctgcaacaattagcaaagcgtctctgcaatattgtgtaaaacaatcagaactaaaggcactttca
aaggctataaatgagccctgtgacacaaactggagatattatggagatagctgctatgggttcttcaggcacaacttaacatgggaagagagtaagcagctgactgacatgaatgctactctctgaagattgacaaccggaacattgtggagtacatcaagccagactcattta
tcgttgggtcggattatctcgcagaagtcgaatgaggtctggaagtgggagatggctcgggtatctcagaaaatattgttgggttttgggaagatgaaaaggaaatatgaattgtgcttattttcataatgggaaatgcacccctacctctgtgagaacaacattattaatgtgtgagagga
aggctggcatgaccaaggtggaccaactaccttaa