

TMEM31-transmembrane protein 31 Gene

Reference: PTXBC029575

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 31; testicular secretory protein Li 58		
For research use only.			

GENE DESCRIPTION	transmembrane protein 31
GENE ID	203562
ORIGIN SPECIES	Human
ACCESSIONS	BC029575
VECTOR	pDONR223
NCBI SYMBOL	TMEM31
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagggttaacagaaaagagtgagggaacaacaactcaagcccaacaactctaatgcaccaatgaagatcaagaagaagaatccaacagtcagaacagcatactccagcaaggcagcgaacacaaagagcagacacacagccatccagatgtcattgccttcacgtaggaca
cctacaacatccagcgacagacgatcaacctcttggaagtcctccgtggcctactgagtggaatttcaacccctatcgattgctctcttttgagctttatctctgaatttctctggtgtttaaagaagcctccatgacatatccattgtctgaaagcccagatggaaagatcggaactgcccc
tcatactccacctcttcgactctccacctctactctacaagttttctctctacaattcttccctttcttttattctctctgtactctctctttttattattgtcttcattctgatctctctga