

C11orf64-chromosome 11 open reading frame 64 Gene

Reference: PTXBC029583

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf64; NCRNA00301; long intergenic non-protein coding RNA 301		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 64
GENE ID	283197
ORIGIN SPECIES	Human
ACCESSIONS	BC029583
VECTOR	pDONR223
NCBI SYMBOL	C11orf64
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccatcatctcacattgtgtgaaagaagaaccagaatgttgggttgatgatgccatggatgaatttaaacagtttggccttgcaataattggtatacttctaattgcatgtgaaattatactcttttaacgaaagatgaaactattcagtgccacatgtaccagtaacagaggagca
aagcaaatcttatactcaagaactgcagctgctgttagatctacctgggtggtccatagagaaactgctcagagaacctgctttacctgcctaa