

C19orf34-chromosome 19 open reading frame 34 Gene

Reference: PTXBC029585

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf34; CSNK1G2 antisense RNA 1		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 34
GENE ID	255193
ORIGIN SPECIES	Human
ACCESSIONS	BC029585
VECTOR	pDONR223
NCBI SYMBOL	C19orf34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggttttagctgcatcaggtgcaaaagagaacacccacaaatcaggagagagggaagccgccacccatggagggaagctgaggagtgaacgccaggatggcaggccggctcctcgagcactgcagacaccagaccaaggcgcactgcatcaccacgagaagg
ccccgtccatggcagaagccaacctcaccgcgcgactgcaccggagacgcccacgaggacagccggcctctgagaacgcagctgagccgtcacaggtcgccatcagcaaaaggcaagagtcagattatgcaaacggaggcaaggctgcatggactccagagcctcagacgcc
atcaaaaaagcaagtgacacatctgcatcaaatccatcacaagacggcctcagactcgggggaaacggagctcgtctgcccacccagactggagtgccctttccctcatcgtgtggaggccgccggcgacatccccgggttccctccactcgaaggcagtcacctc
ggttgccagctcctgctctgtgatgacccacgagggtgcagagtgaaaggcacagcggggttcggcctaactccccggggagccctaggcgagagggttggtcctggggagctgccaggctcttatatccctggcctccatccagtcgccgacgactcagctcctgcccct
gtgcacggggcagcgtctgtag