

TMEM217-transmembrane protein 217 Gene

Reference: PTXBC029657

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf128; dJ355M6.2; transmembrane protein 217		
For research use only.			

GENE DESCRIPTION	transmembrane protein 217
GENE ID	221468
ORIGIN SPECIES	Human
ACCESSIONS	BC029657
VECTOR	pDONR223
NCBI SYMBOL	TMEM217
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgcaggatgttctcctcatgggtggcatctctctcttaataccaccagttcttcatcttgacctgaaccagaagacacacatttctatgaggccaagttcagcatctacgtggactcaaagtcggagctagtcacttggaccctgttccacagggcctaataatcagcactggcct
ctccctcaccaccatcatcgcgtgcttctctcttattgtatccacaagaatatctacatgggctgctgatctatgccatgtggatcatctacagctcatcaacttctccatagtcctgctcctcaacgggatcatcaagatcacttcaagacgctgagttattgcactggatcttccaa
atctcacatgttctgcactttttctgtctgcccctcatcgtcaagcatgcatacaacctttacaaggaatccagactgtggcaggaacgccgccacaggctctgtctccaccttgcagtgaacttatga