

ELF5-E74-like factor 5 (ets domain transcription factor) Gene

Reference: PTXBC029743

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ESE2; ETS-related transcription factor Elf-5; E74-like factor 5 (ets domain transcription factor); epithelium-restricted ESE-1-related Ets factor; epithelium-specific Ets transcription factor 2; E74 like ETS transcription factor 5		

For research use only.

GENE DESCRIPTION	E74-like factor 5 (ets domain transcription factor)
GENE ID	2001
ORIGIN SPECIES	Human
ACCESSIONS	BC029743
VECTOR	pDONR223
NCBI SYMBOL	ELF5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttggactcgggtgacacacagcaccttctgcctaatgcatccttctgcgatccctgatgtcgtgactgatctgttcagcaatgaagagtactacccctgctttgagcatcagacgcctgtgactcactgacatcagtcacccctgaatactggactaagcgccatgtgtgggagtg
gctccagttctctgcgaccagtacaagttggacaccaattgcatctcttctgcaactcaacatcagtgccctgcagctgtgagcatgacacagagaggagttctcgaggcagctggcctctgcggcgagtgactgtacttcctccagaacatccgcacacaaggttactccttttta
atgacgctgaagaagcaagccacatcaagactatgctgattccaactccttgaacaagtgccatcaaaagtcagactgtcacagtcatagtagaacaagcctccaaagttctcatctatgggaattgtacgagacctgcttctatctcctgaagaaactgtggcattctggaatg
ggaagatagggaacaaggaattttcgggtgttaaatcggaagccctggcaaatgtgggacaaggaagaaatgacagaatgacatatgaaagttgagcagagccctgagatactactataaaacaggaaatttggagcgggtgaccgaaggttagtgtacaaatttgaaa
aaatgcacacgggtggcaggaagacaagctatga