

C9orf71-chromosome 9 open reading frame 71 Gene

Reference: PTXBC029780

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein C9orf71; C9orf71; transmembrane protein 252		
For research use only.			

GENE DESCRIPTION	chromosome 9 open reading frame 71
GENE ID	169693
ORIGIN SPECIES	Human
ACCESSIONS	BC029780
VECTOR	pDONR223
NCBI SYMBOL	C9orf71
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagacagaactggcctcattctctgtctcttgcctcctgatgggttctctgatgggtcctggggcctctcattctctgggctccatattcgactgtcaggggagcctgattgcggcctatttgcctctgggggttctgatcctctcagtggaatttctggagcaactatcg
ccaggtgactgaagcaaggagtggtgagggcacatgctccgacaacacctgctcatggggcctgcccgtggccacagtagacaggccagactttaccctccagcttatgaagagagccttgaggtggaaaagcagagctgtcctgcagagagagagcctctggcattctccac
ctctatatacagagacggcctggaattccaggatggaatgactcccaccagaggccccaccatcttatagagagtcctatagccggcctgggtggtgacagcaatctcagaggacgcccagaggcgaaggccaagagtgctga