

FLJ42875-hypothetical LOC440556 Gene

Reference: PTXBC029785

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	long intergenic non-protein coding RNA 982		
For research use only.			

GENE DESCRIPTION	hypothetical LOC440556
GENE ID	440556
ORIGIN SPECIES	Human
ACCESSIONS	BC029785
VECTOR	pDONR223
NCBI SYMBOL	FLJ42875
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggagggcgccggggccggccacccagaggaaaagccagcggagccgggtgggcaggagctggcgccggcggtgcgctggggccttctgacctggatcctccgggtcgccaccgctccctgtcagattcgaaccagattcggtcccaaggctcccacggctc
cagggtgcccagctgtgagccagctcgtcggcggccggatggaggaaagcgtcactggcgctacctctgctgtttctggcaccagagattgagactgtgaacctgaacctcccttgatctcctccggtggggccaggaggagccgaatccaggcccaaccgggctcactgag
aagggtgcaggaaagccccaacgactcctactcgtgtgctggccatactggccaggacgctgcagctggccgagccgctgggtgtctgcaacggggagctccccccgctgggggcaataa