



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

MS4A12-membrane-spanning 4-domains, subfamily A, member 12 Gene

Reference: PTXBC029793

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Ms4a10; membrane-spanning 4-domains subfamily A member 12; membrane spanning 4-domains A12		

For research use only.

GENE DESCRIPTION	membrane-spanning 4-domains, subfamily A, member 12
GENE ID	54860
ORIGIN SPECIES	Human
ACCESSIONS	BC029793
VECTOR	pDONR223
NCBI SYMBOL	MS4A12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatactcatccacgcaacaagccgatctgaagttaaatgaaacacataccacccctaccaccacaagcagctcttatgtgctctgatttcaacagcctctgggttcaatcaacttgaaaaacaaagctcagggtgctcagcgtctcagccctacggcatcatctctcgggaattcttgcta
gcagtcaacccgggtcgaaggaaatatacaaatgataaatccaagtgtgggaacagcagtaataactttaaagaagaagcaaaaggcactagggtgatccagatcatggttgatgatgcacatgtgtttggaattgtttgttttaatatctctcttttagaagaattaggtttgcctcta
ctcgtgttatgtggataccattctgggtgctcttcttattatctcggctctctctctgtcagcatccaaaggagcttccctgtctgtctggaagcagacctgggaatgaacattgttagttctatctggcccttactggagtgatctcgtcgtgtggaatgtgcacaatgggtag
ctggccaagactagctgggcgtgtcttctggaaaaggcattacgccagcgtgatgatctctctctctctgtagctgtgacacagcccaatttggcaaccaagaacaacaccacaactaatatgtctctcgtgttattccaatatgtatgaagcaacccctgtgacaccagcgtct
cttcagctctccctccagatgcacaacactactcagtaatagccctcaataatag

More info: Please contact us for any question: contact@proteogenix.fr