

MGC39900-thymosin beta15b Gene

Reference: PTXBC029803

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Tbeta15b; thymosin beta-15B; thymosin beta 15B		

For research use only.

GENE DESCRIPTION	thymosin beta15b
GENE ID	286527
ORIGIN SPECIES	Human
ACCESSIONS	BC029803
VECTOR	pDONR223
NCBI SYMBOL	MGC39900
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgacaaaccagactgtcggaagtggagaagttgacaggtcaaaactgaagaaaactaatactgaagaaaaatactctccctcaaggaaaaggagtcgtctgtcgccaggctggagtcagcgacacgattcggtcactgcaaactcgcttcccaggtcaagca
gttctcctgcctcagctcccagtagctgggacctcaggcacatcaccaagccaggctaa