

IL34-interleukin 34 Gene

Reference: PTXBC029804

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf77; IL-34; interleukin-34; interleukin 34		
For research use only.			

GENE DESCRIPTION	interleukin 34
GENE ID	146433
ORIGIN SPECIES	Human
ACCESSIONS	BC029804
VECTOR	pDONR223
NCBI SYMBOL	IL34
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccccgggcttcacgtggctgcgtatcttgggatcttctggcgtggccttgggaatgagccttggagatgtggccttgacgcagaatgaggagtgacgtgcacgggttttctgcgggacaagctgcagtacaggagccgacttcagcatgaaacactacttcccatcaac
tacaagatcagtgctgtacgaggggtgttcagaatcgccaacgtcaccaggtgcagagggccaggtgagcagagcgggagctcggtatctgtgggtcttggtagcctcagtgccactgagtcgggtcaggacgtgctcgcaggccaccatcctggaagtacctgcagg
aggtgcagacgtgctgctgaatgtccagcagggcctcacggatgtgaggtcagcccaaggtggaatcgtgtgtccctcttgaatgccccagggccaaacctgaagctgtgcggcccaagccctgctggacaactgctccgggtcatggagctgctgtactgctcctgctga
aacaagctccgtcctaaactggcaggactgtgaggtgccaagtctcagcttgcagccagagccctcattgcagtatgcggccaccagctgtaccctcgcggcgggtgtccccagctccccgctcactccacgggctcgtgagggccggtcagggcacagggcgaggccct
cttgcctga