



Reference: PTXBC029820

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GKP2; GKTA; glycerol kinase 2; ATP:glycerol 3-phosphotransferase 2; glycerokinase 2; glycerol kinase pseudogene 2; glycerol kinase testis specific 2; testis tissue sperm-binding protein Li 77m		

For research use only.

GENE DESCRIPTION	glycerol kinase 2
GENE ID	2712
ORIGIN SPECIES	Human
ACCESSIONS	BC029820
VECTOR	pDONR223
NCBI SYMBOL	GK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]