

C16orf73-chromosome 16 open reading frame 73 Gene

Reference: PTXBC029829

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf73; gs129; meiosis-specific with OB domain-containing protein; meiosis specific with OB domains		

For research use only.

GENE DESCRIPTION	chromosome 16 open reading frame 73
GENE ID	254528
ORIGIN SPECIES	Human
ACCESSIONS	BC029829
VECTOR	pDONR223
NCBI SYMBOL	C16orf73
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagcaactgtaattcctcaaaaccattattacaactaatccagatataccagaagctaacattctgctgaattttatagagaaaataaagaacaaatgttctggatgatgaaattgacagttattcaagaatccataaattaaagtacaatggtgatgtcacacagttgaacaattaaagg
gaaaagctttgaagaatgaaggaagctgaccttctctatggcatctttatgcctacattccacactcaacattgatgatgaaactacaaaagtagttcgaatagatgttcagctgtggttatattgtaaatgaagcatctaacatgtgcacaacttgcaacaaaaactcctggactttaaat
ctgtctttctcagtttccatgtgctgattgatctgactgatcacacaggcacccttcattctgtgctcaggaagtggtgctgaggagactttgggctgcacgttcgttctatcacagagcaaggagtggtgattgaaaattagtgactctctgcaagcttgcagatcctactgaggcaagc
agaaactgtctggacaaaacatgttaa