

MGC34800-hypothetical protein MGC34800 Gene

Reference: PTXBC029861

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	putative uncharacterized protein MGC34800		
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC34800
GENE ID	162137
ORIGIN SPECIES	Human
ACCESSIONS	BC029861
VECTOR	pDONR223
NCBI SYMBOL	MGC34800
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggacctgcccggcgccagagcccaagcggccagcaccgagagctggagacgtcgtgagcaccaggtggcgccggcgcggttcctacaggagcctggagcaggctgcccgggtggctgtgcccggagactggagcggctcctggctccgcaggca
gggccccagagacagctcatgagcagatccgagccccgggccccagcaccgcccggctgcccctcgtcagtcctccgagcgtgccacctggccagccagctgccccctgaacctgcccggcgctccggtctcctgggccccagccaaccccgcg
gctcaggctcccagttcctgtccagcctgcgcccgtccagcggccccggctcctcagctgctctccgactcctatccttggccgccccagcgcggaaccgaccgcgacctaccgccaacatcgcgggtctcgtctcgcagccttttggcgtccgcctcagcgttga