

CARTPT-CART prepropeptide Gene

Reference: PTXBC029882

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CART; cocaine- and amphetamine-regulated transcript protein; CART prepropeptide		
For research use only.			

GENE DESCRIPTION	CART prepropeptide
GENE ID	9607
ORIGIN SPECIES	Human
ACCESSIONS	BC029882
VECTOR	pDONR223
NCBI SYMBOL	CARTPT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagagctcccggtgaggtgctgccccctctggcgccgcccctgctgatgctacctctgttgggtaccggtgccagggaggacgcccagctccagcccgagccctggacatctactctgccgtggatgatgcctccacgagaaggagctgacgaagcgtgcaagaag
tcttgaagaagctcaagagtaaacgtgttcccatctatgagaagaagtatggccaagtcctccatgtgtgacgcccgtgagcagtggtgagtgaggaaaggggcaaggatcgggaagctgtgtgactgtccccgaggaaacctctgcaattcttctctgaagtgttatga