

MS4A5-membrane-spanning 4-domains, subfamily A, member 5 Gene

Reference: PTXBC029884

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD20-L2; CD20L2; TETM4; membrane-spanning 4-domains subfamily A member 5; CD20 antigen-like 2; testis-expressed transmembrane protein 4; testis-expressed transmembrane-4 protein; membrane spanning 4-domains A5		

For research use only.

GENE DESCRIPTION	membrane-spanning 4-domains, subfamily A, member 5
GENE ID	64232
ORIGIN SPECIES	Human
ACCESSIONS	BC029884
VECTOR	pDONR223
NCBI SYMBOL	MS4A5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattcaagcaccgcacacagtcgggtgtttctgtatttctccagaaatcactgttcagaatatgagtcacagaacttccagcagacctttcaactcaaaagcccttgcaaaaattatttctagaaaaatgaaaatcttagggactatccagatcctgtttggaattatgacctttcttt
tggagttatcttctttcacctgttfaaaaccataccaaggttccctttatatttcttccaggaatcattctggggctctgtttgttcattatctggagccttctaattgcagtgaaaagaaaaaacacagaaactctgataatattgagccgaataatgaatttcttagtgccctgggagcaat
agctggaatcattctctcacatttgggtttcatcctagatcaaaactacatttgggttattctcaccaaaatagtcagtgtaaggctgttactgtctgttcttgggaatttgaattacattgatgacttccagcattatgaattatcatttctgcctttctcaatttgggggtgccactcagaggatttga
ttgtgaacaatgtgttga