

RGS5-regulator of G-protein signaling 5 Gene

Reference: PTXBC030059

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MST092; MST106; MST129; MSTP032; MSTP092; MSTP106; MSTP129; regulator of G-protein signaling 5		

For research use only.

GENE DESCRIPTION	regulator of G-protein signaling 5
GENE ID	8490
ORIGIN SPECIES	Human
ACCESSIONS	BC030059
VECTOR	pDONR223
NCBI SYMBOL	RGS5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcaaaggactgcagctttgccccactcatgctgaaaggccaaggagattaagatcaagttgggaattctctccagaagccagactcagttggtgaccttgcattccgtacaatgagaagccagagaaccagccaagaccagaaaacctcgtggacgaggccctgcag
tggcgtgattccctggacaaactctgcagaacaactatggacttccagtttcaaaagtctcgaagtcgaattcagtgaggaaaacctgagttctggtgattcctgtgaggattacaagaagatcaagtcacctgccaagatggctgagaaggcaaaagcaaatatgaagaattcattca
aacggaggctcctaaaggaggtgaattatgaccacttactaaggacatcacatgaagaacctgggtgaaccttccctgagcagcttgacatggccagaaaagaatccatgccctgatggaaaaggattctctgacctctgtgctgctctgagtttatcaggagttatcaagtag