

SH3BGRL3-SH3 domain binding glutamic acid-rich protein like 3 Gene

Reference: PTXBC030135

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SH3BGRL3-like protein; HEL-S-297; SH3BP-1; TIP-B1; SH3 domain-binding glutamic acid-rich-like protein 3; SH3 domain binding glutamic acid-rich protein like 3; SH3 domain-binding protein 1; TNF inhibitory protein; epididymis secretory protein Li 297; SH3 d		

For research use only.

GENE DESCRIPTION	SH3 domain binding glutamic acid-rich protein like 3
GENE ID	83442
ORIGIN SPECIES	Human
ACCESSIONS	BC030135
VECTOR	pDONR223
NCBI SYMBOL	SH3BGRL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcggcctgcgcgtctacagcagctcggcaccggtcccgcaaatcaagtcaccagcagcaggtgacccgaatcctggatgggaagcgcataaccagctagtggacatctccaggacaacgccctgagggatgagatgcgagccttggcaggcaacccaaggc
caccacccccagattgtcaacggggaccagctactgtgggactatgagctctctgtgaggtgtggaacaaacacgctgcaggagtctctgaagctggcttga