



Antibodies and Proteins

## ProteoGenix SAS

## Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

[www.ProteoGenix.science](http://www.ProteoGenix.science)

**GBAP-glucosidase, beta; acid, pseudogene Gene**

Reference: PTXBC030240

|                               |      |                     |                |
|-------------------------------|------|---------------------|----------------|
| Product type:                 | DNA  | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      |      |                     |                |
| <b>For research use only.</b> |      |                     |                |

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | glucosidase, beta; acid, pseudogene                                                               |
| <b>GENE ID</b>          | 2630                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC030240                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | GBAP                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

atggccagcctgctgacttctcatccgcacactacactatgcagacacccctgatgatttcacgttgcaaaccttcagctcccagaggaaagataccaagctcaagatacccttgatccaggagccctcaggtggcccgagctccggtttcaactcttgcagacccctggacatcacccactc  
gtccacacacaaaggagcggggaaatgggaaggccgcatcaaggcagaccccgagacatctacacacagacgtggccagatcatgtgaagttctgactgatctgctgagacacgtttacgttctggcgatgacagctgaataatgagctctgtggttgatgctg  
atccctctccagctcgggctctaccctcgaaatcagcagagactctatgctgctgacttacccttcccaacgggtatccacacaaagtctgctgactgactgacacacgttgcgtgccctgggcaaaagtgtgtactgacgacgcccaagagcagctctaaag  
gcctgtgtgggttccagattctgggagcagatgtgtgcggtacgctctgtggatcggggatgcagtcacagccacatcaaacacctctgtaccatgtgtgctgcgtgacccgactggaacccatctgtagacatcaacaaagcacacttttcaaacagcccatgttaccac  
tctgcccatctcagcaagttctactctgagcgtccggacgagtggtggctggttgcgcagtcagaagaacgacccggacgcagtggtgactgatccatccgatgctctctgtgtgtgctgctctaaacgcctcctctaaaggatgtgccttaccatcaaggaatcctgtgggtctctcg  
agacactctcactctctactcactcactcactcactctgctgcctcccaatga