

C1orf211-chromosome 1 open reading frame 211 Gene

Reference: PTXBC030279

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf211; NCRNA00337; long intergenic non-protein coding RNA 337		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 211
GENE ID	148645
ORIGIN SPECIES	Human
ACCESSIONS	BC030279
VECTOR	pDONR223
NCBI SYMBOL	C1orf211
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaattgaaaagagatttaatttttaaaatcaaaagtcactcagctgcagtcacatgggaaggagcctggatgtgaatgttcctgggagctgagagctcctgtgactgtccagtttctgtggtctccacagactcaggacatctccctctgctccagacacctcaggaccac
cagagggcactcatcgaccttcaccctggaaggagctgtcctgggcccacccactgtgctgggaagccctctggtctctctctcaatgtatga