

KIAA0513-KIAA0513 Gene

Reference: PTXBC030280

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein KIAA0513		
For research use only.			

GENE DESCRIPTION	KIAA0513
GENE ID	9764
ORIGIN SPECIES	Human
ACCESSIONS	BC030280
VECTOR	pDONR223
NCBI SYMBOL	KIAA0513
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagacccagaggtccccgtgggtcgtactgactttggcctgagggcaccacccctctccctggaggcaccacccctgtgctgagggacggcgtgctcctggggacgggtgcatcagagagtgagaccactgagctcgggacagtgagaatgacatgggcgagtcgcccctgcacccgtctctgggaccaagaccgctgtctctccaacgagtccttctctccaacagagcaccgagtcctaccaggatgaagagacccctggcactcagggtctcatgctggtctacgtggagaagatctctctggggggaggacttggatcaggaggagaaagccaagttggagagtgactgagcagtgaaaatggaaaaggccgggagtggtttgctcgtatcgtgagtgcccagcgtgcaactccaagtgtgtctcagaggcaactctaccgctggtgagctctttgagtggtgctgttcgagtgatcatcagatggatgactttggcctgccaagaaacctatgaccatgtcttcaactactaccacatcgaaaaccacagctgtgccccggagtcgccgggagagcccgccggcagcctcgtactcctgaaatccgcaaacagctgggtggccgaaaagaaggacatgccgagcggctgctgaagaacacctggccagagactgagaatgtcaagggtcttctgggggctggagaccaagctgaaggggccccctggccaggaggaacgaggaagacgagaacaacccaggagagaagcggccagggctgtgaccgctacagccccgagggacgaaaagaagggggagaagatctactgtacacgcacctgaa gcaacagcccatctggtaa