

FLJ21438-hypothetical protein FLJ21438 Gene

Reference: PTXBC030281

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAS protein activator like-3; RAS protein activator like 3		

For research use only.

GENE DESCRIPTION	hypothetical protein FLJ21438
GENE ID	64926
ORIGIN SPECIES	Human
ACCESSIONS	BC030281
VECTOR	pDONR223
NCBI SYMBOL	FLJ21438
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgtctccactgccccggccaggtccactccagctctccgaggggagaagccggctctgccccccgggacgtcccaagcacaccctctcatctccaagagccagttctgcgcagcgttcgcgctcagagagttggcccgccacggccggacgaagagcg
gccccctggcgccgccccgggtgcagcgacgcagagtgctccggtccgctgctgccccgcccgaatctgcggggccctggccgcgacccaaaggctccctgagcatgggaccagcggcccgccggccttgaccgggactccgctcgtgctcgggaagc
cgtcggtaccctggcagcgccaaatggaccagccgcaagaccgaaccaggcactgggcacgcaccgacctgtgaacaagtgtggcagagctgcagtgccaggtggccgctctgctgagggagcagaagtgtctgctccgctcgtggagtcgctgagcaccaaatccggcctt
gacggagcagcaggagcagctgcccggccagctgcaggatctgactccagctccgtgctgggagctcagagtttatttcagagcacacctaacaagcaatgaaggccacagctctgaaaaacctggagcaccgctaatgagatggagagaactcaggtcagctgagggatg
ctgtccagagcctgcagcttttccaaggacgcgggggtcttgagtgcaacccagccccctcaagcacctgctcaatggagacaccactga