

STAG3L1-stromal antigen 3-like 1 Gene

Reference: PTXBC030509

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DKFZp434A0131; FLJ78546; MGC131759; MGC40269; STAG3L2; STAG3L3; OTTHUMP00000209277; OTTHUMP00000209427; STAG3-like protein 1; stromal antigen 3-like protein 1; stromal antigen 3-like protein 2; stromal antigen 3-like 1		

For research use only.

GENE DESCRIPTION	stromal antigen 3-like 1
GENE ID	54441
ORIGIN SPECIES	Human
ACCESSIONS	BC030509
VECTOR	pDONR223
NCBI SYMBOL	STAG3L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtggaagagagcaacgccagagccaggtgccagaggactttctcagcttctgtctctctttgtggagagcaagctccacgaccacgtcttacttagtagacaacctgtgggactgtgcagggactcagctgaaggactgggaggtctgacaagcctgctgctggag
aaggaccagagcagctgccatggagccagggccagggaccttcacctcctagggtga