



Reference: PTXBC030532

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARCI11; HAI; MT-SP1; MTSP1; PRSS14; SNC19; TADG15; TMPRSS14; suppressor of tumorigenicity 14 protein; membrane-type serine protease 1; prostamin; serine protease 14; serine protease TADG-15; suppression of tumorigenicity 14 (colon carcinoma); suppression		

GENE DESCRIPTION	suppression of tumorigenicity 14 (colon carcinoma)
GENE ID	6768
ORIGIN SPECIES	Human
ACCESSIONS	BC030532
VECTOR	pDONR223
NCBI SYMBOL	ST14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]